

HOTLINE: Effective May 18, 2020

2002109

Protein Electrophoresis with Reflex to Immunofixation, Serum

SPEP REFLEX

Performed: Sun-Sat
Reported: 1-5 days

Specimen Required: Collect: Serum Separator Tube (SST).

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1.5 mL serum to an ARUP Standard Transport Tube. (Min: 1.5 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma.

Stability (collection to initiation of testing): After separation from cells: Ambient: Unacceptable; Refrigerated: 1 week; Frozen: 1 month