
2002871 PML-RARA Detection by RT-PCR, Quantitative**PML QNT**

Specimen Required: Collect: Whole blood or bone marrow in lavender (EDTA).

Specimen Preparation: **Whole Blood:** Transport 5 mL whole blood. (Min: 3 mL)

Bone Marrow: Transport 3 mL bone marrow. (Min: 1 mL)

Refrigerate immediately. Specimens must be received within 48 hours of collection due to lability of RNA.

Storage/Transport Temperature: **Whole Blood and Bone Marrow:** CRITICAL REFRIGERATED. Separate specimens must be submitted when multiple tests are ordered.

Unacceptable Conditions: Serum, plasma, ambient or frozen bone marrow or whole blood, CSF, or FFPE tissue. Specimens collected in anticoagulants other than EDTA. Severely hemolyzed or clotted specimens.

Ambient whole blood and ambient bone marrow specimens past 7 days will be cancelled. Refrigerated whole blood or bone marrow past 7 days will be canceled.

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 48 hours; Frozen: Unacceptable