

TEST CHANGE

Coxsackie A Antibodies (Serotypes 2, 4, 7, 9, 10 and 16), Serum
2002932, COX A AB

Specimen Requirements:

Patient Preparation:

Collect: Plain red or serum separator tube (SST).

Specimen Preparation: Transfer 2 mL serum to an ARUP standard transport tube. (Min: 1 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Room temperature. Also acceptable: Refrigerated or frozen.

Unacceptable Conditions:

Remarks:

Stability: Ambient: 1 week; Refrigerated: 2 weeks; Frozen: 1 month

Methodology: **Semi-Quantitative** Complement Fixation

Performed: Varies

Reported: **6-11**~~3-8~~ days

Note:

CPT Codes: 86658 x6

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report