

TEST CHANGE

Coxsackie A Antibodies (Serotypes 2, 4, 7, 9, 10, and 16), Serum
2002932, COX A AB

Specimen Requirements:

Patient Preparation:

Collect: Plain red or serum separator tube (SST).

Specimen Preparation: Transfer 2 mL serum to an ARUP standard transport tube. (Min: 1 mL)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: ~~Frozen.~~ ~~Room temperature.~~ Also acceptable: ~~Room temperature~~ ~~Refrigerated~~ or ~~refrigerated~~ ~~frozen.~~

Unacceptable Conditions:

Remarks:

Stability: Ambient: 1 week; Refrigerated: 2 weeks; Frozen: 1 month

Methodology: Semi-Quantitative Complement Fixation

Note:

CPT Codes: ~~86317~~ ~~86658~~ x6

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

~~Refer to~~ ~~By~~ report