
2003128

Tapentadol and Metabolite, Urine, Quantitative

TAPENTA UR

Performed: Mon
Reported: 1-8 days

Specimen Required: Collect: Random urine.

Specimen Preparation: Transfer 2 mL urine with no additives or preservatives to an ARUP Standard Transport Tube. (Min: 1 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Specimens exposed to repeated freeze/thaw cycles.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: 3 years.