

HOTLINE: Effective November 15, 2021

2004247 CEBPA MUT **CEBPA** Mutation Detection

Specimen Required: Collect: Lavender (EDTA) or bone marrow (EDTA).

Specimen Preparation: Whole Blood: Do not freeze. Transport 5 mL whole blood (Min: 1 mL)

Bone Marrow: Transport 3 mL bone marrow. (Min: 1 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma, serum, FFPE tissue blocks/slides, or frozen tissue, DNA extracted by a non-CLIA lab. Specimens

collected in anticoagulants other than EDTA or sodium heparin. Clotted or grossly hemolyzed specimens. Stability (collection to initiation of testing): Ambient: 24 hours; Refrigerated: 5 days; Frozen: Unacceptable