

HOTLINE: Effective May 17, 2021

2005018

Celiac Disease (*HLA-DQ2*, and *HLA-DQ8*) Genotyping

HLA CELIAC

Specimen Required: Collect: Lavender (EDTA), pink (K₂EDTA), or yellow (ACD Solution A or B).

Specimen Preparation: Transport 3 mL whole blood. (Min: 1 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Frozen specimens in glass collection tubes.

Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: 1 month