

TEST CHANGE

Chymotrypsin, Fecal

2005160, CHYMOTRYP

Specimen Requirements:

Patient Preparation: Patients receiving pancreatic enzymes should discontinue taking the enzymes at least 5 days prior to stool collection.

Collect: Random stool.

Specimen Preparation: Transfer 1 g stool to a clean, unpreserved stool transport vial (ARUP Supply #40910). Available online through eSupply using ARUP Connect^(TM) or contact ARUP Client Services at (800-) 522-2787. (Min: 0.5 g) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: ~~Frozen~~Refrigerated. Also acceptable: Refrigerated~~Frozen~~.

Unacceptable Conditions:

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: 2 weeks; Frozen: 1 month

Methodology: Quantitative Enzymatic Assay/Spectrophotometry

Performed: Varies

Reported: ~~8-14~~3-18 days

Note:

CPT Codes: 84311

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report