

---

**2005419****Androstenediol Glucuronide Quantitative, Serum or Plasma****ANDRO GLUC****Performed:** Varies**Reported:** 8-18 days**Specimen Required:** Collect: Plain Red, Lavender (EDTA), or Serum Separator Tube (SST).Specimen Preparation: Separate from cells within **45 minutes**. Transfer 1 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.5 mL)**Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.**Storage/Transport Temperature: Frozen. Also acceptable: Room temperature or refrigerated.Stability (collection to initiation of testing): Ambient: 6 days; Refrigerated: 6 days; Frozen: 3 years