

TEST CHANGE

MPL Mutation Detection by Capillary Electrophoresis

2005545, MPL

Specimen Requirements:

Patient Preparation:

Collect: Whole blood or bone marrow in lavender (EDTA) or pink (K2EDTA). ~~+~~

Specimen Preparation: Whole Blood: Do not freeze. Transport 5 mL whole blood. (Min: 1 mL)
Bone Marrow: Transport 3 mL bone marrow. (Min: 1 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma, serum, FFPE tissue blocks/slides, or frozen tissue.
Specimens collected in anticoagulants other than EDTA.
Clotted or grossly hemolyzed specimens.

Remarks:

Stability: Refrigerated: 7 days; Frozen: Unacceptable

Methodology: Capillary Electrophoresis

Note: The test will detect *MPL* mutations W515K, W515L, W515A, and S505N.

CPT Codes: 81338

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Refer to report.

Reference Interval: