

TEST CHANGE

Natalizumab Antibodies

2005593, NATAL ABS

Specimen Requirements:

Patient Preparation:

Collect: Plain red. Also acceptable: Serum separator tube (SST).

Specimen Preparation: Allow blood to clot at room temperature for 30 minutes. Separate serum from cells within 1 hour. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Frozen. Also acceptable: Refrigerated

Unacceptable Conditions:

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: 2 Weeks; Frozen: 2 Weeks

Methodology: Qualitative Bridging Enzyme-Linked Immunosorbent Assay

Note:

CPT Codes: 83516

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report

HOTLINE NOTE: There is a price change associated with this test. Please contact ARUP Client Services at (800) 522-2787 for additional information.