

TEST CHANGE

Serotonin Release Assay (Heparin Dependent Platelet Antibody), Unfractionated Heparin
2005631, SRA UFH

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube.

Specimen Preparation: Transfer 5 mL serum to ARUP Standard Transport Tubes. (Min: 1 mL)

Transport Temperature: Frozen. Also acceptable: Refrigerated.

Unacceptable Conditions:

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: 1 week; Frozen: **2 years**~~Indefinitely~~.

Methodology: Qualitative Serotonin Release Assay

Note:

CPT Codes: 86022

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

Reference Interval:

<u>Test Number</u>	<u>Components</u>	<u>Reference Interval</u>
	SRA, Unfractionated Heparin	Negative

[Refer to report](#)