

TEST CHANGE

Glutathione Total

2006328, GLUT TOT

Specimen Requirements:

Patient Preparation:

Collect: Yellow (ACD solution B) (ARUP supply #49658) or yellow (ACD solution A). Available online through eSupply using ARUP Connect(TM) or contact ARUP Client Services at 800-522-2787.

Specimen Preparation: Transport whole blood in original collection container. ACD solution B: 10 mL (Min: 8.5) or ACD solution A: 8.5 mL (Min: 6.5 mL) or transfer 1 mL whole blood from a fully drawn tube to an ARUP standard transport tube (Min: 1 mL).
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Critical refrigerated

Unacceptable Conditions:

Remarks: Collection tube must be filled greater than 50% of total volume.

Stability: Ambient: Unacceptable; Refrigerated: 3 weeks; Frozen: Unacceptable

Methodology: Quantitative Kinetic Assay

Note:

CPT Codes: 82978

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

[Refer to report](#) [By Report](#)