

HOTLINE: Effective **October 5, 2020**

2006328

Glutathione Total

GLUT TOT

Specimen Required: Collect: Yellow top (ACD Solution B) (ARUP supply #49658) or Yellow (ACD Solution A). Available online through eSupply using ARUP Connect™ or contact ARUP Client Services at (800) 522-2787.

Specimen Preparation: Transport 10 mL or 8 mL blood in original collection container. (Min: 8.5 mL or 6.5 mL).

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Storage/Transport Temperature: **Critical Refrigerated**

Unacceptable Conditions: Hemolyzed specimens

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 3 weeks; Frozen: Unacceptable