

HOTLINE: Effective July 5, 2022

2006332

Exome Sequencing, Trio

EXOME SEQ

Specimen Required: Collect: Lavender (EDTA) or yellow (ACD Solution A or B). Peripheral blood required. Contact ARUP's genetic counselor at (800) 242-2787 ext. 2141 if there are questions prior to test submission.

AND Maternal Specimen: Lavender (EDTA) or yellow (ACD Solution A or B). Peripheral blood required.

AND Paternal Specimen: Lavender (EDTA) or yellow (ACD Solution A or B). Peripheral blood required.

New York State Clients: ARUP cannot facilitate testing for New York patients. Please work directly with a New York-approved laboratory.

Specimen Preparation: Patient Specimen: Transport 3 mL whole blood. (Min: 1 mL)

AND Maternal Specimen: Transport 3 mL whole blood. (Min: 1 mL)

AND Paternal Specimen: Transport 3 mL whole blood. (Min: 1 mL)

Storage/Transport Temperature: Refrigerated.

Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable