

HOTLINE: Effective November 18, 2019

2006450 Hepatitis Delta Antigen by ELISA**HEPD AG****Performed:** Varies**Reported:** 7-10 days**Specimen Required:** Collect: Plain Red or Serum Separator Tube (SST).Specimen Preparation: Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)Storage/Transport Temperature: Frozen. Also acceptable: Refrigerated.Unacceptable Conditions: Grossly hemolyzed or lipemic specimens.Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 1 week; Frozen: 3 months