
2006450 Hepatitis Delta Antigen by ELISA**HEPD AG**

Specimen Required: Collect: Plain Red or Serum Separator Tube (SST).

Specimen Preparation: Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Storage/Transport Temperature: **CRITICAL FROZEN.**

Unacceptable Conditions: Grossly hemolyzed or lipemic specimens.

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: **Unacceptable**; Frozen: **2 weeks**