
2006550 Thyroglobulin by LC-MS/MS, Serum or Plasma

**THYROG
MS**

Specimen Required: Collect: Serum separator tube or green (sodium or lithium heparin), Potassium EDTA
Specimen Preparation: Separate from cells: Transport 1.5 mL serum or plasma. (Min: 0.7 mL)
Storage/Transport Temperature: Refrigerated or Frozen.
Unacceptable Conditions: Samples left ambient for greater than 1 day; **Grossly lipemic samples.**
Stability (collection to initiation of testing): After separation from cells: Ambient: 1 day; Refrigerated: 1 week; Frozen: 1 year