



HOTLINE: Effective May 18, 2020

---

2007136

von Willebrand Factor Collagen Binding

VWF C BIND

**Specimen Required:** Collect: Light Blue (CTAD).

Specimen Preparation: Transfer 0.5 mL citrated plasma to an ARUP Standard Transport Tube. (Min: 0.3 mL)

Storage/Transport Temperature: **CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.**

Unacceptable Conditions: Hemolyzed specimens.

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 2 weeks

**CPT Code(s):** 83520