

HOTLINE: Effective September 6, 2022

2007136

von Willebrand Factor (VWF) Collagen Binding

VWF C BIND

Performed: Varies
Reported: 3-10 days

Specimen Required: Collect: Light blue (sodium citrate).
Specimen Preparation: Transfer 1 mL citrated plasma to an ARUP Standard Transport Tube. (Min: 0.5 mL)
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.
Storage/Transport Temperature: **CRITICAL FROZEN.**
Unacceptable Conditions: Hemolyzed specimens.
Stability (collection to initiation of testing): Ambient: 24 hours; Refrigerated: 24 hours; Frozen: 3 weeks