

HOTLINE: Effective March 7, 2022

---

|         |                                                                     |            |
|---------|---------------------------------------------------------------------|------------|
| 2007370 | Periodic Fever Syndromes Panel, Sequencing and Deletion/Duplication | PRFEVERPAN |
|---------|---------------------------------------------------------------------|------------|

---

**Methodology:** Massively Parallel Sequencing  
**Performed:** Varies  
**Reported:** 3 weeks

**Specimen Required:** Collect: Lavender or pink (EDTA) or yellow (ACD solution A or B).  
Specimen Preparation: Transport 3 mL whole blood. (Min: 2 mL)  
Storage/Transport Temperature: Refrigerated.  
Unacceptable Conditions: Serum or plasma; grossly hemolyzed or frozen specimens; saliva; buccal brush or swab, FFPE tissue.  
Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable

**HOTLINE NOTE:** Remove information found in the Remarks field.