

HOTLINE: Effective May 2, 2022

---

2007567	Luteinizing Hormone (LH), Pediatric	LH PEDIA
---------	-------------------------------------	----------

---

**Performed:** Varies  
**Reported:** 5-13 days

**Specimen Required:** Collect: Serum separator tube (SST). Also acceptable: Lavender (EDTA).  
Specimen Preparation: Separate serum from cells within 45 minutes. Transfer 1 mL serum to an ARUP Standard Transport Tube.  
(Min: 0.5 mL) Also acceptable: Plasma.  
**Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.**  
Storage/Transport Temperature: Frozen.  
Stability (collection to initiation of testing): Ambient: 24 hours; Refrigerated: 48 hours; Frozen: 6 months