
2007578

High Molecular Weight Kininogen (HMWK), Activity

HIGH MOLE

Specimen Required: Patient Prep: Do not draw from an arm with a heparin lock or heparinized catheter.

Collect: Light Blue (Sodium Citrate).

Specimen Preparation: Transfer 2 mL plasma to an ARUP Standard Transport Tube. (Min: 1 mL)

Storage/Transport Temperature: Frozen.

Stability (collection to initiation of testing): Ambient: 4 hours; Refrigerated: Unacceptable; **Frozen: 2 weeks**