

TEST CHANGE

Kell Antigen Typing - Patient

2007731, K AG

Specimen Requirements:

Patient Preparation:

Collect: Lavender (K2EDTA) or **pPink** (K2EDTA).

Specimen Preparation: Do not freeze. Transport 7 mL whole blood. (Min: 0.5 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Separator tubes.

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: 1 week; Frozen: Unacceptable

Methodology: Qualitative Agglutination~~Hemagglutination~~

Note: When the K (KEL1) antigen is positive, the k (KEL2) antigen may be performed at an additional charge.

CPT Codes: 86905

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Refer to~~By~~ report