

HOTLINE: Effective May 17, 2021

---

2007912

Ethyl Glucuronide Screen with Reflex to Confirmation, Urine

ETG SCR

**Specimen Required:** Collect: Random urine.

Specimen Preparation: Transfer 4 mL urine with no additives or preservatives to an ARUP Standard Transport Tube. (Min: 1 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions:

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: 1 month