

2007914

EPOR Mutation Detection by Sequencing**EPOR****Specimen Required:** Collect: Lavender (EDTA).Specimen Preparation: **Whole Blood:** Do not freeze. Transport 5 mL whole blood. (Min: 1 mL)Storage/Transport Temperature: Refrigerated.Unacceptable Conditions: Plasma, serum, FFPE tissue blocks/slides, or frozen tissue, DNA extracted by a non-CLIA lab, bone marrow. Specimens collected in anticoagulants other than EDTA or sodium heparin. Clotted or grossly hemolyzed specimens.Stability (collection to initiation of testing): Ambient: 24 hours; Refrigerated: 5 days; Frozen: **Unacceptable**