

Quarterly HOTLINE: Effective June 4, 2018

2008326 Hydrocarbon and Oxygenated Volatiles Panel, Blood

HYDRO OX B

Note: Acetaldehyde is an unstable compound post-collection and will both form and degrade under certain sample handling conditions. Even when extreme precautions are taken to maintain the integrity of Acetaldehyde during sample collection, transport and analysis, the results will be affected under typical collection and laboratory procedures.

Test includes: Benzene, Ethylbenzene, Styrene, Toluene, Xylenes (o,m,p), n-Heptane, n-Hexane, Methylpentanes (2- and 3- Isomers), Pentane, n-Butanol, Ethanol, Isopropanol, n-Propanol, Methanol, Acetaldehyde, Acetone, Methyl Ethyl Ketone, Methyl Isobutyl Ketone, Methyl n-Butyl Ketone, Ethyl Acetate, Diethyl Ether, and Methyl Tertiary Butyl Ether.

HOTLINE NOTE: There is a component change associated with this test.

Remove component 0091149, Methyl Acrylate