

TEST CHANGE

NPS Stimulants (Bath Salts ~~Qualitative~~ Panel), Urine

2008650, BATH SALTS

Specimen Requirements:

Patient Preparation:

Collect: Urine.

Specimen Preparation: Transport 12 mL urine in an ARUP standard transport tube. ~~Standard Transport Tube~~. (Min: 0. 47 mL)
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Frozen. Also acceptable: Refrigerated.

Unacceptable Conditions:

Remarks:

Stability: Ambient: 1 week; Refrigerated: 2 weeks~~1 month~~; Frozen: 2 weeks~~1 month~~

Methodology: Semi-Quantitative~~Qualitative~~ High Performance Liquid Chromatography-Tandem Mass Spectrometry

Note:

CPT Codes: 80371 ~~(Alt code: G0480)~~

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report