

HOTLINE: Effective July 6, 2021

2008650 Bath Salts Qualitative Panel, Urine BATH SALTS

Methodology: Qualitative High Performance Liquid Chromatography-Tandem Mass Spectrometry (HPLC-MS/MS)

Specimen Required: Collect: Urine.

Specimen Preparation: Transport 2 mL urine in an ARUP Standard Transport Tube. (Min: 0.7 mL)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Storage/Transport Temperature: Frozen. Also acceptable: Refrigerated.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: 1 month