
2008682 Anabolic Steroids, Urine - Screen with Reflex to Confirmation

STEROIDS

Specimen Required: Collect: Urine.

Specimen Preparation: **Transfer** 4 mL urine to ARUP Standard Transport Tubes. (Min: 1.6 mL)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Storage/Transport Temperature: Refrigerated. Also acceptable: Frozen.

Stability (collection to initiation of testing): Ambient: **Unacceptable**; Refrigerated: 1 month; Frozen: 1 month