

---

**2008893****Amphiphysin Antibody, IgG****AMPHIPHYS****Performed:** Mon, Thu, Sat  
**Reported:** 1-4 days**Specimen Required:** Collect: Serum separator tube.Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.30 mL)Storage/Transport Temperature: Refrigerated.Unacceptable Conditions: Plasma. Contaminated, heat-inactivated, hemolyzed, or lipemic specimens.Stability (collection to initiation of testing): After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 month**CPT Code(s):** 84182