

TEST CHANGE

Antimicrobial Level - Rifabutin by HPLC, Serum or Plasma 2009363, RIFABU

Plain red. Also acceptable: Green (sodium heparin).
Separate from cells ASAP or within one hour of collection. Transfer 2 mL serum or plasma to an ARUP standard transport tube and freeze immediately. (Min: 0.5 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.
CRITICAL FROZEN.
Severely hemolyzed or thawed specimens.
Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 1 month
Quantitative High Performance Liquid Chromatography (HPLC)
Varies
3-10 13-17 days
80299
Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Effective Date: February 3, 2025