

TEST CHANGE

Giardia lamblia Antibodies Panel by ELISA

2009410, GIA PAN

Specimen Requirements:

Patient Preparation:

Collect: Plain red or serum separator tube (SST).

Specimen Preparation: Remove serum from cells within one hour. Transfer 1 mL serum to an ARUP ~~standard transport tube~~ **Standard Transport Tube** and freeze immediately. (Min: 0.5 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Frozen

Unacceptable Conditions:

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: 1 week; Frozen: 5 weeks

Methodology: ~~Semi~~-Quantitative Enzyme-Linked Immunosorbent Assay

Performed: Varies

Reported: ~~6-13~~**3-10** days

Note:

CPT Codes: 86674 x3

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Reference Interval:

By report