

Quarterly HOTLINE: Effective November 12, 2018

2010075 Amphetamines, Urine, Quantitative

AMPS UR

Reference Interval:

Effective November 11, 2018

Drugs Covered	Cutoff Concentrations
Amphetamine	50 ng/mL
Methamphetamine	200 ng/mL
Methylenedioxymphetamine (MDA)	200 ng/mL
Methylenedioxymethamphetamine (Ecstasy, MDMA)	200 ng/mL
Methylenedioxyethylamphetamine (Eve, MDEA)	200 ng/mL
Phentermine	200 ng/mL

Interpretive Data:

Methodology: Quantitative Liquid Chromatography-Tandem Mass Spectrometry.

Positive cutoff: 200 ng/mL unless specified below:

Amphetamine 50 ng/mL

For medical purposes only; not valid for forensic use.

The absence of expected drug(s) and/or drug metabolite(s) may indicate non-compliance, inappropriate timing of specimen collection relative to drug administration, poor drug absorption, diluted/adulterated urine, or limitations of testing. The concentration value must be greater than or equal to the cutoff to be reported as positive. Interpretive questions should be directed to the laboratory.

See Compliance Statement B: www.aruplab.com/CS

CPT Code(s): 80325; 80359 (Alt code: G0480)

HOTLINE NOTE: There is a component change associated with this test.

Add component 2010082, Phentermine, Urn, Quant