
2010673

CALR (Calreticulin) Exon 9 Mutation Analysis by PCR

CALR

Specimen Required: Collect: Whole blood or bone marrow (EDTA).

Specimen Preparation: **Whole Blood:** Do not freeze. Transport 5 mL whole blood. (Min: 1 mL)

Bone Marrow: Do not freeze. Transport 3 mL bone marrow. (Min: 1 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma, serum, FFPE tissue blocks/slides, or frozen tissue. Specimens collected in anticoagulants other than EDTA. Clotted or grossly hemolyzed specimens.

Stability (collection to initiation of testing): Ambient: 24 hours; Refrigerated: 5 days; Frozen: Unacceptable

Extracted DNA: Ambient: 1 month; Refrigerate: Indefinitely; Frozen: Indefinitely