

HOTLINE: Effective July 1, 2019

2010808	Human Immunodeficiency Virus Type 1 (HIV-1) Drug Resistance (PhenoSense GT Plus Integrase)	PHENO PLUS
----------------	---	-------------------

Specimen Required: Collect: Lavender (EDTA) or plasma preparation tube (PPT).

Specimen Preparation: Separate plasma from cells within 6 hours of collection. Transfer 5 mL plasma to an ARUP Standard Transport Tube and freeze immediately. (Min: 3 mL)

Storage/Transport Temperature: **CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.**

Remarks: Provide patient's most recent viral load and viral load collection date.

Unacceptable Conditions: Thawed specimens.

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: Unacceptable; **Frozen: 2 weeks**