

2010905

Collagen Type VII Antibody, IgG by ELISA

COLLAG 7



Immunodermatology Required Clinical
Information Form (Serum)

Methodology: Semi-Quantitative Enzyme-Linked Immunosorbent Assay
Performed: Varies
Reported: 7-14 days

Specimen Required: Collect: Plain red or serum separator tube (SST).
Specimen Preparation: Transfer 2 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)
Storage/Transport Temperature: Refrigerated.
Unacceptable Conditions: Hemolyzed or lipemic specimens. Plasma.
Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 2 weeks; Frozen: Indefinitely

Interpretive Data: Refer to report