

TEST CHANGE

Beryllium Lymphocyte Proliferation, Blood

2011052, BE LYM PRO

Specimen Requirements:

Patient Preparation:

Collect: Green (sodium heparin).

Specimen Preparation: Specimen must be received at performing laboratory with 24 hours of collection. Do not send to ARUP Laboratories. For direct submission instructions, please contact ARUP Referral Testing at (800) 242-2787, ext. 5145. Gently invert several times to mix and prevent clotting. ~~Do~~ do not centrifuge. Transport 40 mL whole blood in the original collection tube(s). (Min: 30 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Send Monday-Thursday only. Specimen must be sent directly to performing laboratory. ~~Room at room~~ temperature.

Unacceptable Conditions: Centrifuged specimens.

Remarks:

Stability: Ambient: 30 hours; Refrigerated: Unacceptable; Frozen: Unacceptable

Methodology: Cell Culture

Performed: Varies

Reported: ~~14-17~~ 15-20 days

Note:

CPT Codes: 86353

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By Report



*A nonprofit enterprise of the University of Utah
and its Department of Pathology*

Effective Date: **February 3, 2025**