

2011172

Urogenital Ureaplasma and Mycoplasma Species by PCR

UR MYCOPCR

Performed: Mon, Wed, Fri
Reported: 2-5 days

Specimen Required: Collect: Genital swab, rectal swab, or urine. Also acceptable: upper respiratory swabs, bronchoalveolar lavage, sputum, and tracheal aspirates.
Specimen Preparation: Transfer genital swab, rectal swab, respiratory swab, or 1 mL urine to viral transport media (ARUP supply #12884) available online through eSupply using ARUP Connect™ or contact ARUP Client Services at (800) 522-2787.
BAL, sputum, or tracheal aspirate: Transfer 1 mL to an empty sterile container. (Min: 0.5 mL)
Storage/Transport Temperature: Frozen
Remarks: Specimen source required.
Stability (collection to initiation of testing): Ambient: 48 hours; Refrigerated: 10 days; Frozen: 14 days