

HOTLINE: Effective February 18, 2020

2011375 Occupation Screen - MMR/VZV Antibody Assessment Panel, IgG

MMRV PAN

Reference Interval:

Components	Reference Interval	
Measles Virus (Rubeola) Antibody IgG	13.4 AU/mL or less	Negative - No significant level of detectable measles (rubeola) IgG antibody.
	13.5-16.4 AU/mL	Equivocal - Repeat testing in 10-14 days may be helpful.
	16.5 AU/mL or greater	Positive - IgG antibody to measles (rubeola) detected, which may indicate a current or past exposure/immunization to measles (rubeola).
Mumps Virus Antibody IgG	8.9 AU/mL or less	Negative - No significant level of detectable IgG mumps virus antibody.
	9.0-10.9 AU/mL	Equivocal - Repeat testing in 10-14 days may be helpful.
	11.0 AU/mL or greater	Positive - IgG antibody to mumps virus detected, which may indicate a current or past exposure/immunization to mumps virus.
Rubella Virus Antibody IgG	Less than 9 IU/mL	Not Detected.
	9-9.9 IU/mL	Indeterminate - Repeat testing in 10-14 days may be helpful.
	10 IU/mL or greater	Detected.
Varizella-zoster Virus Ab IgG	Effective February 18, 2020	
	134.9 IV or less	Negative - No significant level of detectable varicella-zoster IgG antibody.
	135.0-164.9 IV	Equivocal - Repeat testing in 10-14 days may be helpful.
	165.0 IV or greater	Positive - IgG antibody to varicella-zoster detected, which may indicate a current or past varicella-zoster infection.

HOTLINE NOTE: There is a numeric map change associated with this test.

Change the numeric map for component 2011399, Varizella-zoster Virus Ab IgG from XXXXX to XXXX.X.