

TEST CHANGE

NPS Stimulants (Bath Salts ~~Qualitative~~ Panel), Serum ~~or~~ Plasma

2011411, BATHSLT SP

Specimen Requirements:

Patient Preparation:

Collect: Plain ~~red, lavender~~ Red, Lavender (K2EDTA), or ~~p~~ Pink (K2EDTA).

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer ~~12~~ mL serum or plasma to an ARUP ~~standard transport tube~~ Standard Transport Tube. (Min: 0.47 mL)
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Frozen.

Unacceptable Conditions: Gel ~~separator tubes~~ Separator Tubes. Thawed specimens.

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: ~~24~~48 hours; Frozen: 2 weeks

Methodology: ~~Semi-Quantitative~~ Qualitative High Performance Liquid Chromatography-Tandem Mass Spectrometry

Note:

CPT Codes: 80371 (~~Alt code: G0480~~)

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report