

TEST CHANGE

Allergen, Food, Cashew IgG

2011727, CASHEW IGG

Specimen Requirements:

Patient Preparation:

Collect: Plain red or serumSerum separator tube (SST).

Specimen Preparation: Separate <u>serum from cells ASAP or within 2 hours of collection.</u>

Transfer 0.5 mL serum or plasma to an ARUP standard transport tube. Standard Transport Tube. (Min: 0.252 mL).

Refer to "Allergen Specimen Collection Instructions" at www.aruplab.com/testing/resources/specimen.)

Effective Date: July 21, 2025

Transport Temperature: Refrigerated.

Unacceptable Conditions: Hemolyzed, icteric, or lipemic specimens.

Remarks:

Stability: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year

Methodology: Quantitative ImmunoCAP Fluorescent Enzyme Immunoassay

Performed: Sun

Reported: 1-8 days

Note:

CPT Codes: 86001

New York DOH Approval Status: Specimens from New York clients will be sent out to a New

York DOH approved laboratory, if possible.

Interpretive Data:

Values less than 2.00 mcg/mL represent absent or undetectable levels of allergen-specific IgG antibody.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

Reference Interval:

Less than 7.19 mcg/mL



Effective Date: July 21, 2025