

TEST CHANGE

Hepatitis E Virus (HEV) Antibodies, IgG and IgM

2012023, HEV PAN

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube (SST). Also acceptable: lavender (K2EDTA K2-EDTA), lavender (K3EDTA K3-EDTA), or pink (K2EDTA).

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 0.5 mL serum or plasma to an ARUP standard transport tube. Standard Transport Tube: (Min: 0.3 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Specimens containing particulate material.

Remarks:

Stability: After separation from cells: Ambient: 48 hours Unacceptable; Refrigerated: 14 days 2 weeks; Frozen: 30 days Indefinitely (avoid repeated freeze/thaw cycles)

Methodology: Qualitative Enzyme-Linked Immunosorbent Assay (ELISA)

Note:

CPT Codes: 86790 x2

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

~~This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA-certified laboratory and is intended for clinical purposes.~~

Reference Interval:

Test Number	Components	Reference Interval
	Hepatitis E Virus Ab, IgG	Negative
	Hepatitis E Virus Ab, IgM	Negative

Deleted Cells

Deleted Cells