
2012039 Lysozyme, Serum**LYSO SER****Specimen Required:** Collect: Serum separator tube (SST).Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.4 mL)Storage/Transport Temperature: Refrigerated.Unacceptable Conditions: Hemolyzed, lipemic, icteric, or contaminated specimens.Stability (collection to initiation of testing): After separation from cells: Ambient: **Unacceptable**; Refrigerated: **5 days**; Frozen: **1 month**.**Reference Interval:**

Effective November 18, 2019

Less than or equal to 2.75 µg/mL