
2012155

**Charcot-Marie-Tooth (CMT) and Related Hereditary Neuropathies, *PMP22*
Deletion/Duplication with Reflex to Sequencing Panel**

CMT REFLEX

Specimen Required: Collect: Lavender (EDTA), or yellow (ACD Solution A or B).

Specimen Preparation: Transport 3 mL whole blood. (Min: 2 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Serum or plasma; grossly hemolyzed or frozen specimens

Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable