



HOTLINE: Effective May 17, 2021

---

2012211

Barbiturates, Urine Screen with Reflex to Quantitation

BARB RFX U

**Specimen Required:** Collect: Random urine.

Specimen Preparation: Transfer 4.0 mL urine with no additives or preservatives to an ARUP Standard Transport Tube. (Min: 2.0 mL)

Storage/Transport Temperature: Room temperature.

Unacceptable Conditions: Breast milk. Pharmaceutical preparation.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: 1 month