TEST CHANGE

Gabapentin, Urine	
2012227, GABAP U	
Specimen Requirements:	
Patient Preparation:	
Collect:	Random urine.
Specimen Preparation:	Transfer 1 mL urine to an ARUP Standard Transport Tube. (Min: 0.6 mL)
Transport Temperature:	Refrigerated.
Unacceptable Conditions:	
Remarks:	
Stability:	Ambient: 1 month; Refrigerated: 1 month; Frozen: 1 month
Methodology:	Quantitative Liquid Chromatography-Tandem Mass Spectrometry
Performed:	Mon, Wed, Sat
Reported:	1- <u>7</u> € days
Note:	
CPT Codes:	80355 (Alt code: G0480)
New York DOH Approval Status:	This test is New York DOH approved.
Interpretive Data:	
Positive cutoff: 5.0 μ g/mL	

For medical purposes only; not valid for forensic use.

The absence of expected drug(s) and/or drug metabolite(s) may indicate non-compliance, inappropriate timing of specimen collection relative to drug administration, poor drug absorption, diluted/adulterated urine, or limitations of testing. The concentration value must be greater than or equal to the cutoff to be reported as a quantitative result. Interpretive questions should be directed to the laboratory.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

Reference Interval:

By report

