

HOTLINE: Effective May 17, 2021

---

2012245

Methadone, Urine Screen with Reflex to Quantitation

METH RFX U

**Specimen Required:** Specimen Preparation: Transfer 4 mL urine with no additives or preservatives to an ARUP Standard Transport Tube. (Min: 2 mL)  
Storage/Transport Temperature: Room temperature.  
Unacceptable Conditions: Breast Milk. Specimens exposed to repeated freeze/thaw cycles.  
Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: 1 month