

HOTLINE: Effective May 17, 2021

2012269

Propoxyphene, Urine Screen with Reflex to Quantitation

PPXY RFX U

Specimen Required: Collect: Random urine.

Specimen Preparation: Transfer 4 mL urine with no additives or preservatives to an ARUP Standard Transport Tube. (Min: 2 mL)

Storage/Transport Temperature: Room temperature.

Unacceptable Conditions: Breast milk.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: **1 month**

Note: If the specimen screens positive, then Confirmation/Quantitation by LC-MS/MS (ARUP test code **3003726**) will be added to confirm result. Additional charges apply.

HOTLINE NOTE: There is a reflexive pattern change associated with this test.

Add reflex to 3003726, Propoxyphene and Metabolite, Urine

Remove reflex from 2010468, Propoxyphene and Metabolite, Urine, Quantitative