

---

2012273

**Buprenorphine, Urine Screen with Reflex to Quantitation**

**BUPR RFX U**

**Specimen Required:** Collect: Random urine.

Specimen Preparation: Transfer 4 mL urine with no additives or preservatives to an ARUP Standard Transport Tube (Min: 2 mL)

Storage/Transport Temperature: Room temperature.

Unacceptable Conditions: Unknown fluids, pharmaceutical preparations, and breast milk. Specimens exposed to repeated freeze/thaw cycles.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: **1 month**